



ISOCS / LabSOCS

**Calibration software for
Gamma Spectroscopy**

Setting up a detector for measurement

► Setup Hardware


◆ Peak Shaping Parameters

- Rise Time
- Flat Top

◆ Pole-zero

◆ Number of Channels

◆ Signal Gain





► Setup Software

◆ Energy Calibration

◆ Shape Calibration

◆ Efficiency Calibration





◆ Select Sources

- Energy Range
- Half-life
- Cascade effects
- Activity

◆ Select Geometries

- Beakers
- Vials

◆ Pay \$\$\$ for these sources



Source-based Efficiency Calibrations



► Cost of Source-based Calibration


- ◆ Purchase of Sources
- ◆ Replacement of Sources
- ◆ Disposal of Sources
- ◆ Licensing for Sources
- ◆ Calibration Program
- ◆ Sample Preparation

Time consuming

► What if you can't make the sample match the calibration

- ◆ Not enough material
- ◆ Unusual material composition
- ◆ Locked in sealed container


Calibration Sample



?

=

Measured Sample




Corrections required


Gamma-ray Interactions in Material

▶ Gamma-ray Interactions are very well understood.

- ◆ Photoelectric absorption
- ◆ Compton Scattering
- ◆ Pair Production




Energy dependence of gamma ray interaction in Germanium



- ▶ We can exploit this knowledge to create efficiency responses based on the physical parameters of the geometry.
- ▶ Use mathematical models to accurately compute the transport of gamma-rays through these geometries.

Field (In-Situ) measurements



Modeling of Laboratory Samples





► How to handle varieties of samples and containers?



ISOCS/LabSOCS: A Generalized Efficiency Computation Method

- ▶ ISOCS/LabSOCS method is designed to quickly and accurately compute efficiencies for a wide range of geometries.
- ▶ Relies on a factory characterization of the intrinsic efficiency response of the detector.
- ▶ Does not require that the customer have a large inventory of calibration sources (although nominal sources for quality control tracking is highly recommended).




Regulatory Acceptance of Modeling Approaches

- ▶ **“Calibration of Germanium Detectors for In-Situ Gamma-ray Measurements”, N42.28-2002 American National Standards Institute, Inc., 1430 Broadway, New York 10018.**
“One such application [of Monte Carlo Methods] is the calculation of the efficiency or response function for an HPGe detector. Using this approach, detectors can be calibrated for a variety of applications using models and simulations.”
- ▶ **“Measuring, Evaluating, and Reporting Radioactive Material in Liquid and Gaseous Effluents and Solid Waste”, U.S. Nuclear Regulatory Commission Guide 1.21 rev. 2 (June 2009):**
“The use of NIST-traceable sources combined with mathematical efficiency calibrations may be applied to instrumentation used for radiochemical analysis (e.g., gamma spectroscopy systems) if employing a method provided by the instrument manufacturer.”
- ▶ **“A Good Practice Guide for the use of Modelling Codes in Non Destructive Assay of Nuclear Materials”, ESARDA Bulletin No. 42 (November 2009) 26.**
- ▶ **[Proposed revision of N42.14-1999. In committee]**
“Calibration and Use of Germanium Spectrometers for the Measurement of Gamma-Ray Emission Rates of Radionuclides”, N42.14-201x American National Standards Institute, Inc., 1430 Broadway, New York 10018:
“The following approaches may be considered for the calibration of the detector efficiency:
 - a) *Measurement of a standardization coefficient for a specific gamma ray and radionuclide by direct comparison with a standard source of known activity;*
 - b) *Measurement of the full-energy peak efficiency as a function of energy;*
 - c) *Calculation of the peak efficiency as a function of energy with the use of Monte Carlo or other calculation techniques.”*


Flexibility of Modeling

- ▶ With modeling, it is possible to rapidly produce geometries that represent many usual shapes for which source standards may not be readily available.




Include Features Beyond just the Sample

- ▶ Complex objects can be created including absorber materials such as spacers, collimators, and positioning rings.





A Word of Caution

- ▶ While mathematical modeling is a significant time and cost saver compared to source-based calibrations, care must still be taken to model the geometry to an appropriately accurate level.
- ▶ Not all regions are equally important.
- ▶ Some critical parameters are:
 - ◆ Distance from Sample to Detector
 - ◆ Attenuation of gamma-rays by intervening materials




Sensitivity of Distance

- ▶ The efficiency for sources close to the detector are very sensitive to the position of the source.
- ▶ As a general rule of thumb: 1 mm \approx 5% change in efficiency for close geometries.
- ▶ Many beakers have complicated bases and an “effective distance” may not be clear.
- ▶ Also note: Reference source standards in epoxy matrices may deform beaker (from Rxn heat) thus changing the geometry.
- ▶ Direct modeling of beaker base provides accurate description of the container.



3D Interactive Geometry Composer

► **Powerful 3D visualization provides immediate feedback for faster geometry development**



Edit dimensions - Complex Pipe

Description: Lead-lined Drum

Comment:

Units: mm cm m in ft

No.	Description	d.1	d.2	d.3	d.4	d.5	Material	Density	Rel. Conc.
1	Pipe	1.125	10.5	17.25	-17.25		iron	7.86	
2	Source 1	10.499	0	17.25	16.125		iron	7.86	0.00
3	Source 2	10.499	0	-16.125	-17.25		iron	7.86	0.00
4	Source 3	1.499	9	16.124	-16.124		lead	11.4	0.00
5	Source 4	8.999	0	16.124	14.75		lead	11.4	0.00
6	Source 5	8.999	0	-14.75	-16.124		lead	11.4	0.00
7	Source 6	8.999	0	5	-14.74		dirt1	1.6	1.00
8	Source 7	0	0	0	0			0	0.00
9	Source 8	0	0	0	0			0	0.00
10	Source 9	0	0	0	0			0	0.00
11	Source 10	3	1	-5	4	30	sand	1.7	4.00
12	Absorber 1	0						0	
13	Absorber 2	0						0	
14	Source - Detector	12	0	0	0	0			

Buttons: OK, Cancel, Apply, View Drawing...

Info
Name: Spherical Source

Styles
Visible:
Transparent:
Color: RGB (0, 0, 255)

Dimensions
Units: cm
d2.1: 15
d2.2: 60
d2.3: 20
d2.4: 20

Material
Identifier: csteel
Density: 7.86
Concentration: 0.00

d2.1 Source diameter

Ready




◆ Detectors, “basic” Collimators, and Templates are rendered in the 3D virtual space.

◆ Data entry is fundamentally similar to the present geometry composer (although not identical).


Templates - sample shape

- ▶ Many standard sample geometries are available for LabSOCS
 - ◆ Box, Cylinder, Sphere, Beaker, Marinelli beaker...
- ▶ Beaker editor for complex shape (cylindrically symmetric)
- ▶ Many templates are also available for ISOCS
 - ◆ Drum, Container, Pipe, Box, surface measurement...
 - ◆ With collimator, housing, or without

Simplified Marinelli Beaker



Material editor




- ▶ **Easy to create:**
 - ◆ New element
 - ◆ List of currently used elements


LabSOCS Beaker Editor

- ▶ Used to create the *.bkr files for the complex beaker template in LabSOCS
- ▶ Simplifies creating of the complex beaker geometries by providing a visual output
- ▶ Created geometry may consist of boundaries for different materials, shields and collimators with a variety of shapes

Beaker View Window



Beaker Edit Window



Geometry report

Geometry Composer Report 


Date: Thursday, May 05, 2011 - 17:26:09
Description: SG15 à 5c0 + support complet
Comment:
File Name: C:\Genie2k\CAMFILES\Cattenom12\CAT2 SG15 det1.geo
Software: LabSOCS
Template: SIMPLIFIED_BEAKER, Version: Custom Beaker (sg15catsupport2.BKR)
Detector: b11084
Environment: Temperature = 22 °C, Pressure = 760 mm Hg, Relative Humidity = 30%
Integrator: Convergence = 1.00%, MDRPN = 2¹ (16), CRPN = 2¹ (16)
Sample Weight: 16.3 g

Dimensions (mm)									
No.	Description	d.1	d.2	d.3	d.4	d.5	d.6	Material	Density
1	Beaker								
2	Sample	40.672						epoxy	1.1
3	Absorber 1								
4	Absorber 2								
5	Source - Detector	49							

List of energies for efficiency curve generation

40.0	59.5	89.0	122.1	165.9	391.7	661.7	898.0
1173.2	1332.5	1836.0					

Custom Simplified Beaker



The diagram shows a 3D perspective view of a beaker with a blue liquid inside. Dimensions are labeled: 2.1 for the height of the liquid, 5.1 for the diameter of the liquid surface, 3.1 for the diameter of the beaker neck, and 4.1 for the diameter of the beaker body. Below the beaker is a yellow rectangular absorber labeled 'Absorber(s)'. To the right is a 2D cross-sectional view of the beaker with a vertical dashed line representing the 'Reference Plane'.

Page 1 of 2

► ISOCS / LabSOCS report

◆ All geometry parameters


- Dimensions
- Materials

◆ Draw

- Template used
- Final result

ISOCS - LabSOCS



► Create an efficiency curve to calculate activities





► Create a geometry file

- ◆ This file can be used to create a new geometry
 - Dimensions modification
 - Materials modification
- ◆ Can be used to cascade summing correction

Modeling for Sample Optimization




- ▶ Cylindrical sample with 500 cc of 1.6 g/cc Soil.
- ▶ Ran multiple efficiency computations with varying container diameter and height but with fixed sample volume.
- ▶ Determine maximum efficiency.
- ▶ All detectors have maximum efficiency when sample diameter is about 12 – 14 cm, regardless of energy.
- ▶ In all cases a high aspect ratio (e.g. BEGe's) detector has greater efficiency than a low aspect ratio detector of similar relative efficiency.



How Certain is your Measurement?

- ▶ If you don't know all the parameters of your measurement, it is possible to use modeling techniques to estimate the uncertainty due to these "not well known" parameters.
- ▶ Canberra provides an ISOCS interface called the ISOCS Uncertainty Estimator (IUE) to estimate these uncertainties
- ▶ The IUE software
 - ◆ Helps the user determine which parameters to concentrate his effort in accurately determining efficiency
 - ◆ A structured and defensible method to quickly create the uncertainty of the efficiency calibration, and to propagate those errors to the final result
 - ◆ A useful investigative tool to evaluate, optimize and choose between various counting choices
- ▶ Different analysis modes:
 - ◆ Sensitivity Analysis: Determine which geometrical parameters are most important to know well (focus effort on minimization of uncertainties of important variables)
 - ◆ Uncertainty Analysis: Put in all known geometrical uncertainties and analyzed to determine the Total Geometric Uncertainties (more reliable accountancy)



▶ Source-Based Efficiency

- ◆ Expensive – Purchase sources/licensing
- ◆ Limited – Few Number of Geometries
- ◆ Labor Intensive – Make Samples Match the Calibration

▶ Mathematical Efficiency

- ◆ Take Advantage of the Well-known Properties of Gamma-ray Physics
- ◆ Generally Applicable
 - In-Situ
 - Non-destructive Assay
 - Laboratory
- ◆ Defendable – Results are reproducible
- ◆ Complex corrections -- Cascade Summing
- ◆ Uncertainty Modeling
- ◆ Measurement Optimization
- ◆ Focus on Making the Calibration Match the Sample

▶ **Mathematical Efficiency modeling not only saves time, but it gives a deeper understanding of the measurement and greater confidence in the results.**